



Child Care Bed Medical Bed for Children or Pediatric Children

Welcome to wholesale and buy customized Child Care Bed Medical Bed for Children or Pediatric Children from our factory. You can rest assured to buy discount Child Care Bed Medical Bed from us. We look forward to cooperating with you, if you want to know more, you can consult us now, we will reply to you in time!

Product Description

Manual Crank Flat Child Care Bed Medical Bed For Children ,Pediatric Children Hospital Beds For Sale

Our crew through skilled training. Skilled expert knowledge, solid sense of service, to meet the service needs of customers for Wholesale price China Popular Fashion Professional China Child Care Bed Medical Bed for Children, Pediatric Children Manufacturers & Suppliers, Your inquiry is going to be highly welcomed and also a win-win affluent growth are what we are anticipating.

Wholesale price China Popular Fashion Professional China Child Care Bed Medical Bed for Children, Pediatric Children, Baby Bed, Strict quality control is executed in each link of the whole production process. We sincerely hope to establish the friendly and mutual-beneficial cooperation with you. Based on high quality products and solutions and perfect pre-sales /after-sales service is our idea, some clients had cooperated with us for more than 5 years.

Features

- 1) Overall size (L*W*H): 2080*950*500mm
- 2) Frame: steel, treated by sand-blasting and electrostatic powder-coating
- 3) Bedboard: steel, 4-part
- 4) Headboard/footboard: ABS, Detachable
- 5) Side rail: Detachable, ABS
- 6) Handset controller: remote, movable
- 7) Motor: Timotion, Linak, L&K
- 8) Castors: four Central-controlled silent wheels, $\phi 125$ mm
- 9) Weight load: 300kg
- 10) Standard accessories: I.V pole 1pc, IV pole holes 4pcs, Drainage hooks 2pcs

Product Display



Packaging&Details

Packaging	20GP: 40PCS 40HQ: 80PCS
Warranty	2 Years
Delivery time	20 Days after Receiving 30% T/T In Advance

Product Details

Model	SAE-YG-01
Overall Size	2080*950*500mm
Bed Frame	Treated By Electro-Coating And Powder-Coating
Headboard / Footboard	ABS head&foot board
Bedboards	Cold-Rolled Steel Bedboards
Handrails	Aluminum Alloy protective railing
Hand Controller	Crank
Bed Base	Steel Frame
Wheels	4 Double - Side Silcent wheels (150mm)
Load Bearing	Fully Tested Sturdy Construction Capable Of Taking A Maximum User Weight Of Up To 300kg

Function

Backrest Max Upward Angle	75°
---------------------------	-----

**FAQ**

1. which products you product?

Medical bed /medical trolley/home care bed / Electric Adjustable Beds

2. How do you ensure the products quality?

Production/inspection report before your order packed. we have obtained ISO Certificate.

3. How about your after-sale service?

We have a professional team to responsible for your after-sale service. Once the order confirmed, we have professional service follow this order.

4. What is your minium order quantity? If only buy one set to check, whether you can offer the best discount price for me?

Our Minium order quantity is one set. we will offer the best competitive price for each customer.

What is your delivery time? if the delivery time is urgent, how do you deal with?

Within 15 working days after order confirmed. If your delivery time is very urgent, we can delivery goods in 7 days. We take customer satisfying as our first goal.

5. What's the payment term?

Usually we accept T/T, L/C, Paypal, and Wester Union.

Warranty: 1-2 years due to different products, and lifetime maintaining.

A. Offer door to door service;

B. We send you free parts by courier and technical guidance by video/call;

C. We support you free training.