



Child Care Bed Medical Bed

Welcome to wholesale and buy customized Child Care Bed Medical Bed from our factory. You can rest assured to buy discount product from us. We look forward to cooperating with you, if you want to know more, you can consult us now, we will reply to you in time!

Product Description

Manual Crank Flat Child Care Bed Medical Bed For Children ,Pediatric Children Hospital Beds For Sale

As a way to ideal meet up with client's desires, all of our operations are strictly performed in line with our motto "High Top quality, Competitive Cost, Fast Service" for Professional Design China Child Care Bed Medical Bed Manufacturers & Suppliers, We are willing to give you the best suggestions on the designs of your orders in a professional way if you need. In the meantime, we keep on developing new technologies and creating new designs so as to make you ahead in the line of this business.

Professional Design China Child Care Bed Medical Bed, Nursing Bed, At present our sales network is growing continually, improving service quality to meet customer's demand. If you are interested in any products and solutions , be sure to contact us at anytime. We are looking forward to forming successful business relationships with you in near future.

Features

- 1) Overall size (L*W*H): 1650*800*500mm
- 2) Frame: steel, treated by sand-blasting and electrostatic powder-coating
- 3) Bedboard: steel, 4-part
- 4) Headboard/footboard: ABS, Detachable
- 5) Side rail: Detachable, ABS
- 6) Handset controller: remote, movable
- 7) Motor: Timotion, Linak, L&K
- 8) Castors: four Central-controlled silent wheels, $\phi 125$ mm
- 9) Weight load: 300kg
- 10) Standard accessories: I.V pole 1pc, IV pole holes 4pcs, Drainage hooks 2pcs

Product Display



Packaging&Details

Packaging	20GP: 40PCS 40HQ: 80PCS
Warranty	2 Years
Delivery time	20 Days after Receiving 30% T/T In Advance

Product Details

Model	SAE-YG-02
Overall Size	1650*800*500mm
Bed Frame	Treated By Electro-Coating And Powder-Coating
Headboard / Footboard	ABS head&foot board
Bedboards	Cold-Rolled Steel Bedboards
Handrails	Aluminum Alloy protective railing
Hand Controller	Crank
Bed Base	Steel Frame
Wheels	4 Double - Side Silcent wheels (150mm)
Load Bearing	Fully Tested Sturdy Construction Capable Of Taking A Maximum User Weight Of Up To 300kg

Function

Backrest Max Upward Angle	75°
Footrest Max Upward Angle	40°

FAQ

1. which products you product?

Medical bed /medical trolley/home care bed / Electric Adjustable Beds

2. How do you ensure the products quality?

Production/inspection report before your order packed. we have obtained ISO Certificate.

3. How about your after-sale service?

We have a professional team to responsible for your after-sale service. Once the order confirmed, we have professional service follow this order.



4. What is your minimum order quantity? If only buy one set to check, whether you can offer the best discount price for me?

Our minimum order quantity is one set. We will offer the best competitive price for each customer. What is your delivery time? If the delivery time is urgent, how do you deal with it? Within 15 working days after order confirmed. If your delivery time is very urgent, we can deliver goods in 7 days. We take customer satisfaction as our first goal.

5. What's the payment term?

Usually we accept T/T, L/C, Paypal, and Western Union.

Warranty: 1-2 years due to different products, and lifetime maintenance.

A. Offer door-to-door service;

B. We send you free parts by courier and technical guidance by video/call;

C. We support you free training.