



Electric Hospital Bed with Five-Function(CPR)

Welcome to wholesale and buy customized Electric Hospital Bed with Five-Function(CPR) from our factory. You can rest assured to buy discount Electric Hospital Bed with Five-Function(CPR) from us. We look forward to cooperating with you, if you want to know more, you can consult us now, we will reply to you in time!

Product Description

Electric Hospital Bed with Five-Function(CPR) made in China can be bought with low price from shuaner. It is a professional high quality Medical Bed manufacturers and factory in China. We are a China company that can provide customized products service. If you want pricelist and quotation, you can ask us by leaving a message. We will get back to you as soon as possible.

Hospital Furniture Five Functions Adjustable Electric ICU Medical Bed

We depend on sturdy technical force and continually create sophisticated technologies to meet the demand of Hot-selling China Good Quality Electric Hospital Bed with Five-Function(CPR), Since the manufacturing unit founded, we have now committed into the progress of new goods. Along with the social and economic pace, we'll continue to carry forward the spirit of "high excellent, efficiency, innovation, integrity", and stay with the operating principle of "credit initially, customer 1st, good quality excellent". We'll produce a excellent foreseeable future in hair output with our companions.

Hot-selling China Good Quality Electric Hospital Bed with Five-Function(CPR), ICU Bed, As an experienced manufacturer we also accept customized order and we could make it the same as your picture or sample specification. The main goal of our company is to live a satisfactory memory to all the customers, and establish a long term business relationship with buyers and users all over the world.

Feature

- StayPut Retractable bed frame design(Patient will feel very convenient while be on drip or used the heart monitor since it won't need to change patient position when the bed back section rise and low.)
- Made of superior cold-rolled steel
- One step forming sleeping bed
- Treated by shot blasting, electrophoresis and powder coating
- Four pieces of PP complete side rails, full protection with zero gap, tuck-away

Product Display



Packaging&Details

Packaging	20GP: 20 PCS 40HQ: 45 PCS
Warranty	2 Years
Delivery time	20 Days after Receiving 30% T/T In Advance

Product Details

Model	A-5b
Overall Size	2190*990*450-730mm
Bed Frame	Treated By Electro-Coating And Powder-Coating
Headboard / Footboard	PP Head / Footboard
Bedboards	Cold-Rolled Steel Bedboards
Handrails	PP Handrails
Hand Controller	Remote
Motor	Quiet And Robust Rlectric Actuators Provide Reliable Operation
	Remote Control Of all Electrically Operated Functions And Movements
	Linak Motor From Denmark And Jiechang From China
Bed Base	Steel Frame
Wheels	4 Double - Side Silcent wheels (150mm)With central locking system
Load Bearing	Fully Tested Sturdy Construction Capable Of Taking A Maximum User Weight Of



	Up To 300kg
--	-------------

Function

Backrest Max Upward Angle	75°
Footrest Max Upward Angle	40°
Height Adjustment	450 - 730 mm
Trendelenburg	15°
Anti-Trendelenburg	15°

Optional Accessories

I.V pole, Nurse control, CPR, Battery, Mattress

FAQ

1. which products you product?

Medical bed /medical trolley/home care bed / Electric Adjustable Beds

2. How do you ensure the products quality?

Production/inspection report before your order packed. we have obtained ISO Certificate.

3. How about your after-sale service?

We have a professional team to responsible for your after-sale service. Once the order confirmed, we have professional service follow this order.

4. What is your minium order quantity? If only buy one set to check, whether you can offer the best discount price for me?

Our Minium order quantity is one set. we will offer the best competitive price for each customer.

What is your delivery time? if the delivery time is urgent, how do you deal with?

Within 15 working days after order confirmed. If your delivery time is very urgent, we can delivery goods in 7 days. We take customer satisfying as our first goal.

5. What's the payment term?

Usually we accept T/T, L/C, Paypal, and Wester Union.

Warranty: 1-2 years due to different products, and lifetime maintaining.

A. Offer door to door service;

B. We send you free parts by courier and technical guidance by video/call;

C. We support you free training.