



Electric Hospital Patient Bed with Five-Function

Shuaner is one of the famous China Electric Hospital Patient Bed with Five-Function manufacturers and suppliers. Our factory specializes in manufacturing of Electric Hospital Patient Bed with Five-Function.

Product Description

Hospital Furniture Five Functions Adjustable Electric ICU Medical Bed

Our enterprise since its inception, usually regards product top quality as business life, repeatedly enhance manufacturing technology, make improvements to product excellent and continuously strengthen enterprise total high quality administration, in strict accordance with all the national standard ISO 9001:2000 for Fixed Competitive Price China Low Price Electric Hospital Patient Bed with Five-Function Manufacturers & Suppliers, President of our company, with whole staff members, welcomes all prospective buyers to go to our corporation and inspect. Let us cooperate hand in hand to create a fantastic long term.

Fixed Competitive Price China Low Price Electric Hospital Patient Bed with Five-Function, Hospital Furniture, With the best technological support, now we have tailored our website for the best user experience and kept in mind your ease of shopping. we ensure that the best reaches you at your doorstep, in the shortest possible time and with the help of our efficient logistical partners i.e DHL and UPS. We promise quality, living by the motto of promising only what we can deliver.

Feature

- StayPut Retractable bed frame design(Patient will feel very convenient while be on drip or used the heart monitor since it won't need to change patient position when the bed back section rise and low.)
- Made of superior cold-rolled steel
- One step forming sleeping bed
- Treated by shot blasting, electrophoresis and powder coating
- Four pieces of PP complete side rails, full protection with zero gap, tuck-away

Product Display



Packaging&Details

Packaging	20GP: 20 PCS 40HQ: 45 PCS
Warranty	2 Years
Delivery time	20 Days after Receiving 30% T/T In Advance

Product Details

Model	Da-5a
Overall Size	2190*990*450-730mm
Bed Frame	Treated By Electro-Coating And Powder-Coating
Headboard / Footboard	PP Head / Footboard
Bedboards	Cold-Rolled Steel Bedboards
Handrails	PP Handrails
Hand Controller	Remote
Motor	Quiet And Robust Rlectric Actuators Provide Reliable Operation
	Remote Control Of all Electrically Operated Functions And Movements
	Linak Motor From Denmark And Jiechang From China
Bed Base	Steel Frame
Wheels	4 Double - Side Silcent wheels (150mm)With central locking system
Load Bearing	Fully Tested Sturdy Construction Capable Of Taking A Maximum User Weight Of Up To 300kg

Function

Backrest Max Upward Angle	75°
Footrest Max Upward Angle	40°
Height Adjustment	450 - 730 mm
Trendelenburg	15°



Anti-Trendelenburg

15°

Optional Accessories

I.V pole, Nurse control, CPR, Battery, Mattress

FAQ

1. which products you product?

Medical bed /medical trolley/home care bed / Electric Adjustable Beds

2. How do you ensure the products quality?

Production/inspection report before your order packed. we have obtained ISO Certificate.

3. How about your after-sale service?

We have a professional team to responsible for your after-sale service. Once the order confirmed, we have professional service follow this order.

4. What is your minium order quantity? If only buy one set to check, whether you can offer the best discount price for me?

Our Minium order quantity is one set. we will offer the best competitive price for each customer.

What is your delivery time? if the delivery time is urgent, how do you deal with?

Within 15 working days after order confirmed. If your delivery time is very urgent, we can delivery goods in 7 days. We take customer satisfying as our first goal.

5. What's the payment term?

Usually we accept T/T, L/C, Paypal, and Wester Union.

Warranty: 1-2 years due to different products, and lifetime maintaining.

A. Offer door to door service;

B. We send you free parts by courier and technical guidance by video/call;

C. We support you free training.