



Electric Medical Bed with Five-Function

Shuaner is one of the famous China Electric Medical Bed with Five-Function manufacturers and suppliers. Our factory specializes in manufacturing of Electric Medical Bed with Five-Function.

Product Description

Electric Medical Bed with Five-Function made in China can be bought with low price from shuaner. It is a professional high quality Medical Bed manufacturers and factory in China. We are a China company that can provide customized products service. If you want pricelist and quotation, you can ask us by leaving a message. We will get back to you as soon as possible.

Hospital Furniture Five Functions Adjustable Electric ICU Medical Bed

The really abundant projects management experiences and 1 to just one provider model make the high importance of business enterprise communication and our easy understanding of your expectations for Hot sale Factory China Electric Medical Bed with Five-Function Manufacturers & Suppliers, Our company has been devoting that "customer first" and committed to helping customers expand their business, so that they become the Big Boss ! Hot sale Factory China Electric Medical Bed with Five-Function, Electric Hospital Bed, We can meet the various needs of customers at home and abroad. We welcome new and old customers to come to consult & negotiate with us. Your satisfaction is our motivation! Allow us to work together to write a brilliant new chapter!

Feature

- StayPut Retractable bed frame design(Patient will feel very convenient while be on drip or used the heart monitor since it won't need to change patient position when the bed back section rise and low.)
- Made of superior cold-rolled steel
- One step forming sleeping bed
- Treated by shot blasting, electrophoresis and powder coating
- Four pieces of PP complete side rails, full protection with zero gap, tuck-away

Product Display



Ningbo Shuaner Medical Equipment Co., Ltd.

www.nbshuaner.com



Packaging&Details

Packaging	20GP: 20 PCS 40HQ: 45 PCS
Warranty	2 Years
Delivery time	30 Days after Receiving 30% T/T In Advance

Product Details

Model	B-5a
Overall Size	2190*980*400-730mm
Bed Frame	Treated By Electro-Coating And Powder-Coating
Headboard / Footboard	PP Head / Footboard
Bedboards	Cold-Rolled Steel Bedboards
Handrails	PP Handrails
Hand Controller	Remote
Motor	Quiet And Robust Rlectric Actuators Provide Reliable Operation
	Remote Control Of all Electrically Operated Functions And Movements
	Linak Motor From Denmark And Jiechang From China
Bed Base	Steel Frame
Wheels	4 Double - Side Silcent wheels (150mm)With central locking system
Load Bearing	Fully Tested Sturdy Construction Capable Of Taking A Maximum User Weight Of Up To 300kg

Function

Backrest Max Upward Angle	75°
Footrest Max Upward Angle	40°
Height Adjustment	400 - 730 mm

Phone: +86-15267896677

E-mail: info@nbshuaner.com



Trendelenburg	15°
Anti-Trendelenburg	15°

Optional Accessories

I.V pole, Nurse control, CPR, Battery, Mattress

FAQ

1. which products you product?

Medical bed /medical trolley/home care bed / Electric Adjustable Beds

2. How do you ensure the products quality?

Production/inspection report before your order packed. we have obtained ISO Certificate.

3. How about your after-sale service?

We have a professional team to responsible for your after-sale service. Once the order confirmed, we have professional service follow this order.

4. What is your minium order quantity? If only buy one set to check, whether you can offer the best discount price for me?

Our Minium order quantity is one set. we will offer the best competitive price for each customer.

What is your delivery time? if the delivery time is urgent, how do you deal with?

Within 15 working days after order confirmed. If your delivery time is very urgent, we can delivery goods in 7 days. We take customer satisfying as our first goal.

5. What's the payment term?

Usually we accept T/T, L/C, Paypal, and Wester Union.

Warranty: 1-2 years due to different products, and lifetime maintaining.

A. Offer door to door service;

B. We send you free parts by courier and technical guidance by video/call;

C. We support you free training.