



Electric Patient Bed with Three-Function

Welcome to wholesale and buy customized Electric Patient Bed with Three-Function from our factory. You can rest assured to buy discount product from us. We look forward to cooperating with you, if you want to know more, you can consult us now, we will reply to you in time!

Product Description

Electric Patient Bed with Three-Function Made in China can be bought with Low Price from Shuaner. It is a professional high quality Products Manufacturers and Factory in China. We can provide Customized Products service. If you want Pricelist and Quotation, you can ask us by leaving a message. We will get back to you as soon as possible.

High Grade 4pcs ABS Handrail ICU Household Nursing Used Three Functions Electric Bed For Hospital sale

1. PP headboard&foot board, 5pcs bars aluminum guiderail.
2. 4pcs 24V DCElectro-motor, 1pc control box.1pc handset
3. 4pcs urine hooks, 1pc IV pole,4pcs luxury castor (2pcs with brake ,2pcs without brake)
4. bed is special design for old patient caring, electric control it like Chair

Our growth depends around the superior machines, exceptional talents and consistently strengthened technology forces for Wholesale Discount China Electric Patient Bed with Three-Function Manufacturers & Suppliers, Our expert complex team will likely be wholeheartedly at your services. We sincerely welcome you to definitely go to our website and company and deliver us your inquiry.

Wholesale Discount China Electric Patient Bed with Three-Function, Hospital Bed, Our company offers the full range from pre-sales to after-sales service, from product development to audit the use of maintenance, based on strong technical strength, superior product performance, reasonable prices and perfect service, we are going to continue to develop, to supply the high-quality merchandise and services, and promote lasting cooperation with our customers, common development and create a better future.

Product Display



Ningbo Shuaner Medical Equipment Co., Ltd.

www.nbshuaner.com



Packaging &Details

Packaging	20GP: 30 PCS 40HQ: 50 PCS
Warranty	2 Years
Delivery time	20 Days after Receiving 30% T/T In Advance

Product details

Model	B-3a
Overall Size	2190*980*400-730mm
Bed Frame	Treated By Electro-Coating And Powder-Coating
Headboard / Footboard	PP Head / Footboard
Bedboards	Cold-Rolled Steel Bedboards
Handrails	PP Handrails
Hand Controller	Remote
Motor	Quiet And Robust Rlectric Actuators Provide Reliable Operation
	Remote Control Of all Electrically Operated Functions And Movements
	Linak Motor From Denmark And Jiechang From China
Bed Base	Steel Frame
Wheels	4 Double - Side Silcent wheels (150mm)With central locking system
Load Bearing	Fully Tested Sturdy Construction Capable Of Taking A Maximum User Weight Of Up To 300kg

Function

Backrest Max Upward Angle	75°
Footrest Max Upward Angle	40°

Phone: +86-15267896677

E-mail: info@nbshuaner.com



Height Adjustment	400- 730 mm
-------------------	-------------

Optional Accessories

I.V pole, Nurse control, CPR, Battery, Mattress

FAQ

1. which products you product?

Medical bed /medical trolley/home care bed / Electric Adjustable Beds

2. How do you ensure the products quality?

Production/inspection report before your order packed. we have obtained ISO Certificate.

3. How about your after-sale service?

We have a professional team to responsible for your after-sale service. Once the order confirmed, we have professional service follow this order.

4. What is your minium order quantity? If only buy one set to check, whether you can offer the best discount price for me?

Our Minium order quantity is one set. we will offer the best competitive price for each customer.

What is your delivery time? if the delivery time is urgent, how do you deal with?

Within 15 working days after order confirmed. If your delivery time is very urgent, we can delivery goods in 7 days. We take customer satisfying as our first goal.

5. What's the payment term?

Usually we accept T/T, L/C, Paypal, and Wester Union.

Warranty: 1-2 years due to different products, and lifetime maintaining.

A. Offer door to door service;

B. We send you free parts by courier and technical guidance by video/call;

C. We support you free training.