



### **Five-Function Electric ICU Hospital Bed**

Welcome to wholesale and buy customized Five-Function Electric ICU Hospital Bed from our factory. You can rest assured to buy discount product from us. We look forward to cooperating with you, if you want to know more, you can consult us now, we will reply to you in time!

#### **Product Description**

Five-Function Electric ICU Hospital Bed made in China can be bought with low price from shuaner. It is a professional high quality Medical Bed manufacturers and factory in China. We are a China company that can provide customized products service. If you want pricelist and quotation, you can ask us by leaving a message. We will get back to you as soon as possible.

Five function Hot Sales hospita Electric Bed for ICU Bed medical furniture  
We constantly think and practice corresponding to the change of circumstance, and grow up.  
We aim at the achievement of a richer mind and body along with the living for Supply  
OEM/ODM China Five-Function Electric ICU Hospital Bed, Currently, we are seeking ahead to even bigger cooperation with abroad buyers determined by mutual benefits. Please sense absolutely free to contact us for more details.

Supply OEM/ODM China Five-Function Electric ICU Hospital Bed, Medical Hospital Bed, We have a dedicated and aggressive sales team, and many branches, catering to our customers. We're looking for long-term business partnerships, and ensure our suppliers that they will undoubtedly benefit in both short and long run.

#### **Product Display**



### Packaging&Details

Packaging	20GP: 35 PCS 40HQ: 65 PCS
Warranty	2 Years
Delivery time	20 Days after Receiving 30% T/T In Advance

### Product Details

Model	E-5b
Overall Size	2160*1080*400-730mm
Bed Frame	Treated By Electro-Coating And Powder-Coating
Headboard / Footboard	PP head&foot board
Bedboards	Cold-Rolled Steel Bedboards
Handrails	PP Handrails
Hand Controller	Remote
Motor	Quiet And Robust Rlectric Actuators Provide Reliable Operation
	Remote Control Of all Electrically Operated Functions And Movements
	Linak Motor From Denmark And Jiechang From China
Bed Base	Steel Frame
Wheels	4 Double - Side Silcent wheels (150mm)
Load Bearing	Fully Tested Sturdy Construction Capable Of Taking A Maximum User Weight Of Up To 300kg

### Function

Backrest Max Upward Angle	75°
Footrest Max Upward Angle	40°
Height Adjustment	400 - 730 mm
Trendelenburg	15°
Anti-Trendelenburg	15°



### **Optional Accessories**

I.V pole, Nurse control, CPR, Battery, Mattress

### **FAQ**

1. which products you product?

Medical bed /medical trolley/home care bed / Electric Adjustable Beds

2. How do you ensure the products quality?

Production/inspection report before your order packed. we have obtained ISO Certificate.

3. How about your after-sale service?

We have a professional team to responsible for your after-sale service. Once the order confirmed, we have professional service follow this order.

4. What is your minium order quantity? If only buy one set to check, whether you can offer the best discount price for me?

Our Minium order quantity is one set. we will offer the best competitive price for each customer.

What is your delivery time? if the delivery time is urgent, how do you deal with?

Within 15 working days after order confirmed. If your delivery time is very urgent, we can delivery goods in 7 days. We take customer satisfying as our first goal.

5. What's the payment term?

Usually we accept T/T, L/C, Paypal, and Wester Union.

Warranty: 1-2 years due to different products, and lifetime maintaining.

A. Offer door to door service;

B. We send you free parts by courier and technical guidance by video/call;

C. We support you free training.