



Five-Function MDF Electric Home Care Bed for Paralysis Patient

Welcome to wholesale and buy customized Five-Function MDF Electric Home Care Bed for Paralysis Patient from our factory. You can rest assured to buy discount product from us. We look forward to cooperating with you, if you want to know more, you can consult us now, we will reply to you in time!

Product Description

5 function electric home care bed and hospital bed for paralysis patient

To be a result of ours specialty and service consciousness, our company has won a superb reputation between customers all around the environment for Rapid Delivery for China Five-Function MDF Electric Home Care Bed for Paralysis Patient, We've been extremely aware of excellent, and have the certification ISO/TS16949:2009. We are dedicated to supply you high quality solutions with acceptable price tag.

Rapid Delivery for China Five-Function MDF Electric Home Care Bed for Paralysis Patient, Electric Medical Home Care Bed, Now, we are trying to enter new markets where we do not have a presence and developing the markets we have the already penetrated. On account of superior quality and competitive price, we will be the market leader, please don't hesitate to contact us by phone or email, if you are interested in any of our products.

Features

- 1) Overall Length: 2150cm
- 2) Overall Width: 1030cm
- 3) Height adjustment: 400-800cm
- 4) Motor: Limoss Motor from German(4 actuators,1 control box,1 remote handset), 2 independent motors on head and foot part,which can adjust the angle of bed height,remote control has a function locking
- 5) Bed surface and bed frame: Cold-Rolled Steel Tube&sheet&plate,thickening,scientific welding,it will not become distorted,make it more stable
- 6) Top coat: Epoxy powder coated finish
- 7) Painting color can be Beige/White or any color
- 8) Head and foot board: anti-fire MDF board,Oak grain panel stucked on the board surface,Fire proof composite board outside,beautiful design,comfortable,home feeling

Product Display



Packaging&Details

Packaging	20GP: 55PCS 40HQ: 108PCS
Warranty	2 Years
Delivery time	20 Days after Receiving 30% T/T In Advance

Product Details

Model	K-5a
Overall Size	2150*1030*400-800mm
Bed Frame	Treated By Electro-Coating And Powder-Coating
Headboard / Footboard	MDF head&foot board
Bedboards	Cold-Rolled Steel Bedboards
Handrails	Aluminum Alloy protective railing
Hand Controller	Remote
Bed Base	Steel Frame
Wheels	4 Double - Side Silcent wheels (150mm)
Load Bearing	Fully Tested Sturdy Construction Capable Of Taking A Maximum User Weight Of Up To 300kg

Function

Backrest Max Upward Angle	75°
Footrest Max Upward Angle	40°
Height Adjustment	400 - 800 mm

FAQ

1. which products you product?

Medical bed /medical trolley/home care bed / Electric Adjustable Beds



2. How do you ensure the products quality?

Production/inspection report before your order packed. we have obtained ISO Certificate.

3. How about your after-sale service?

We have a professional team to responsible for your after-sale service. Once the order confirmed, we have professional service follow this order.

4. What is your minium order quantity? If only buy one set to check, whether you can offer the best discount price for me?

Our Minium order quantity is one set. we will offer the best competitive price for each customer.

What is your delivery time? if the delivery time is urgent, how do you deal with?

Within 15 working days after order confirmed. If your delivery time is very urgent, we can delivery goods in 7 days. We take customer satisfying as our first goal.

5. What's the payment term?

Usually we accept T/T, L/C, Paypal, and Wester Union.

Warranty: 1-2 years due to different products, and lifetime maintaining.

A. Offer door to door service;

B. We send you free parts by courier and technical guidance by video/call;

C. We support you free training.