



High Quality Aged People Easily Controlled Lightweight Electric Power Wheelchair

Welcome to wholesale and buy customized High Quality Aged People Easily Controlled Lightweight Electric Power Wheelchair from our factory. You can rest assured to buy discount Controlled Lightweight Electric Power Wheelchair from us. We look forward to cooperating with you, if you want to know more, you can consult us now, we will reply to you in time!

Product Description

High Quality Aged People Easily Controlled Lightweight Electric Power Wheel chair made in China can be bought with low price from shuaner. It is a professional high quality Medical Bed manufacturers and factory in China. We are a China company that can provide customized products service. If you want pricelist and quotation, you can ask us by leaving a message. We will get back to you as soon as possible.

(Shuaner N-40B)Max Load aluminum framefolding portable electric power wheelchair We attempt for excellence, provider the customers", hopes to be the most beneficial cooperation team and dominator enterprise for staff, suppliers and shoppers, realizes value share and continuous advertising for Reasonable price for China High Quality Aged People Easily Controlled Lightweight Electric Power Wheelchair Manufacturers & Suppliers, Presently, we are looking forward to even larger cooperation with overseas clients based on mutual benefits. Make sure you come to feel no cost to make contact with us for more information. Reasonable price for China High Quality Aged People Easily Controlled Lightweight Electric Power Wheelchair, Electric Wheelchair, We win many reliable customers by rich experience, advanced equipments, skilled teams, strict quality control and best service. We can guarantee all our merchandise. Customers' benefit and satisfaction are always our biggest goal. Be sure to contact us. Give us a chance, give you a surprise.

Product Display



Product Details

Frame material	Aluminum alloy frame
Seat back cushion	Double washable 3D anilox
top speed	1-8km programmable
Charging time	5-12 hours
Climbing angle	≤12 degrees
Rear wheel dimensions	12 inch pneumatic tire (PU solid tire optional), aluminum alloy wheel hub
Front wheel dimensions	7-inch Pu solid tire
Handrail	Rear lift
brake	Electromagnetic brake
Turning radius	830mm
drive	Brush motor 250W * 2 (brushless 250W * 2350W * 2 optional)
Battery	Lithium battery 6 / 12 / 20Ah
Endurance	8/15/25 KM
controller	360 degrees of freedom
Seat height	470mm
Seat depth	420mm
Seat width	460mm
Net weight	16.5kg
Gross weight	18kg
bearing	≤100kg
Unfold size	1000*590*920mm
Folding size	730*330*730mm



Carton Size	750*350*750mm
Optional Accessories	Headrest available
colour	silver

FAQ

1. which products you product?

Medical bed /medical trolley/home care bed / Electric Adjustable Beds

2. How do you ensure the products quality?

Production/inspection report before your order packed. we have obtained ISO Certificate.

3. How about your after-sale service?

We have a professional team to responsible for your after-sale service. Once the order confirmed, we have professional service follow this order.

4. What is your minium order quantity? If only buy one set to check, whether you can offer the best discount price for me?

Our Minium order quantity is one set. we will offer the best competitive price for each customer.

What is your delivery time? if the delivery time is urgent, how do you deal with?

Within 15 working days after order confirmed. If your delivery time is very urgent, we can delivery goods in 7 days. We take customer satisfying as our first goal.

5. What's the payment term?

Usually we accept T/T, L/C, Paypal, and Wester Union.

Warranty: 1-2 years due to different products, and lifetime maintaining.

A. Offer door to door service;

B. We send you free parts by courier and technical guidance by video/call;

C. We support you free training.