



Medical Hydraulic Emergency Transfer Folding Stretcher

Welcome to wholesale and buy customized Medical Hydraulic Emergency Transfer Folding Stretcher from our factory. You can rest assured to buy discount product from us. We look forward to cooperating with you, if you want to know more, you can consult us now, we will reply to you in time!

Product Description

SAE-TC-03 Medical hydraulic emergency transfer folding stretcher

We insist over the principle of improvement of 'High good quality, Effectiveness, Sincerity and Down-to-earth working approach' to deliver you with excellent company of processing for Hot-selling China Ce ISO Medical Hydraulic Emergency Transfer Folding Stretcher Manufacturers & Suppliers, Our corporation has been devoting that "customer first" and committed to helping shoppers expand their small business, so that they become the Big Boss !

Hot-selling China Ce ISO Medical Hydraulic Emergency Transfer Folding Stretcher, Medical Equipment, Profession, Devoting are always fundamental to our mission. We've got always been in line with serving customers, creating value management objectives and adhering to the sincerity, dedication, persistent management idea.

Product Display



Accessories

With pushing handle on the head end, With a crank to operate the foot section, With gas spring for adjust back section

Side rails: 2pcs, Aluminum alloy, collapsible

Back section is operated by gas spring



With the middle castor for steer direction, Bumper castors on 4 corners

Detail information

Functions

Back section up and down :0-75°(by air spring)

Leg section up and down :0-45°(by manual crank)

Height adjust:63-93cm (hydraulic pump)

Trendelenburg and anti-trendelenburg:0-10°(hydraulic pump)

Central brake

Back section available for X-Ray

Advantages

Two seperated Hydraulic pump

-It is adopted two hydraulic pumps;Two kinds of pumps can be chosen:Chinese pumps, Imports pumps.

it is used for the height adjust and the trendelenburg and anti-trendelenburg of the stretcher surface.

Central brake system

-The flat carrier adopt the interantional advanced central-brake system ,it is stable and reliable .

Also ,

the universal /directional castors are configured ,so , one person can be easy operated

Advanced design

-From the overall design and manufacture ,the flat carrier can satisfy your requirement with more humanization,

-it is ideal choice in hospital for carrying patients .

-The main frame structure of flat carrier adopts the high-quality mild steel powder coated.

-it is beautiful and firm for its appearance .

Available for X-Ray

-the back section is made of special electric panel which is imported from USA, so it is available for X-Ray.

Accessories

1 pair of aluminum alloy rails

4 pcs of 5" castors

2 pcs of Hydraulic pump

1 pcs of table

2 pcs of corner bumbers.

FAQ

1. which products you product?

Medical bed /medical trolley/home care bed / Electric Adjustable Beds



2. How do you ensure the products quality?

Production/inspection report before your order packed. we have obtained ISO Certificate.

3. How about your after-sale service?

We have a professional team to responsible for your after-sale service. Once the order confirmed, we have professional service follow this order.

4. What is your minium order quantity? If only buy one set to check, whether you can offer the best discount price for me?

Our Minium order quantity is one set. we will offer the best competitive price for each customer.

What is your delivery time? if the delivery time is urgent, how do you deal with?

Within 15 working days after order confirmed. If your delivery time is very urgent, we can delivery goods in 7 days. We take customer satisfying as our first goal.

5. What's the payment term?

Usually we accept T/T, L/C, Paypal, and Wester Union.

Warranty: 1-2 years due to different products, and lifetime maintaining.

A. Offer door to door service;

B. We send you free parts by courier and technical guidance by video/call;

C. We support you free training.