



New Style Patient Transport Stretcher with Manual Crank Patient Transport Trolley

Welcome to wholesale and buy customized New Style Patient Transport Stretcher with Manual Crank Patient Transport Trolley from our factory. You can rest assured to buy discount product from us. We look forward to cooperating with you, if you want to know more, you can consult us now, we will reply to you in time!

Product Description

New Style Patient Transport Stretcher with Manual Crank Patient Transport Trolley made in China can be bought with low price from shuaner. It is a professional high quality Medical Bed manufacturers and factory in China. We are a China company that can provide customized products service. If you want pricelist and quotation, you can ask us by leaving a message. We will get back to you as soon as possible.

New style patient transport stretcher with manual crank patient transport trolley
Top quality comes first; services is foremost; organization is cooperation" is our small business philosophy which is regularly observed and pursued by our company for Good User Reputation for China New Style Patient Transport Stretcher with Manual Crank Patient Transport Trolley, Together with our efforts, our goods have won the trust of buyers and been very salable equally here and abroad.

Good User Reputation for China New Style Patient Transport Stretcher with Manual Crank Patient Transport Trolley, Transport Stretcher, Our products are sold to the Middle East, Southeast Asia, Africa, Europe, America and other regions, and are favorably appraised by clients. To benefit from our strong OEM/ODM capabilities and considerate services, please contact us today. We will sincerely create and share success with all clients.

Model Number	SAE-TC-02
Overall Size	1950*700*500-900mm
Back Rest	0-75°
Weight Capacity	250KGS
Height Adjustment	630-930mm
Frame	PP Plastic material and painted steel
Platform	Bed board made of PP material

Product Display



Features

Product Size: 1950*700*500-900mmmm

Material: The bed is made of galvanized steel tube and surface coated by electrostatic spray. The back board is made of X-ray revealed American electrotyped board with x-ray cassette. It is equipped with the imported LINAK motor as rise-and-fall system.

Description: The international advanced centrally controlled brake system is applied to the cart and $\Phi 200$ dustproof castors make the cart move steadily, reliably and lightly. The bed can be raised and lowered, tilted forward and backward which is controlled by the treadles lying on two sides of the bed. The bottom cover is made of ABS material and shaped by once injection that can be easily cleaned. The cart surface is of separate structure and the upper part is supported by the imported large hydraulic pole to make the bed never sway. It is completed with infusion stand, drainage hook and collapsible aluminum alloy side rails. It is the best choice for emergency room, ICU, ICU bed and emergency bed for transferring patient.

Technical Features

Full Length 212CM

Full Width 80CM

Bed Base Height Range Maximum Height 93CM

Minimum Height 63CM

Bed Base Angle Position Adjustable

Back section panel 0-75°

Leg section panel 0-45°

Trendelenburg 0-15°



Anti-Trendelenburg 0-15°

Maximum safe load 250KG

Diameter of Centrally Control Silent and Dustproof Castor 20CN

Special Allocation

- Linak Electrical Motor
- Electronic Control Panel
- I.V.Pole
- Rise-and-Fall Guide Wheel
- Mattress (optional)
- Aluminum Alloy Side Rail
- Back Rest may X-ray (Electric Panel Made in America)
- CPR
- Top Coat: Electrostatic Power Coating
- Bed Frame: Aluminum Alloy Side Rail

Using Instruction and Notes

1. Mainly made of PP or ABS plastic and steel-plastic, aluminum and stainless steel, easy to clean and corrosion-proof.
2. PP side rail , bed board , aluminum alloy connecting rails, guide wheel and central control castors.
3. I .V pole , central control mute castros.
4. Hydraulic stretcher cart for operation room transporting

Material: The cart is made of imported PP material. The imported hydraulic rise-and-fall system is used for it.

Description: The international advanced centrally controlled brake system is applied to the cart and $\phi 200$ mm dustproof castors make the cart move steadily, reliably and lightly. The bed can be raised and lowered, tilted forward and backward which is controlled by the treadle lying on two sides of the bed. The bottom cover is made of PP material and shaped by once injection that can be easily cleaned.

The cart surface is of separate structure and the upper part is supported by the imported large hydraulic pole to make the bed never sway. It is completed with infusion stand, drainage hook and collapsible aluminum alloy side rails. It is the best choice for emergency room, ICU, ICU bed and emergency bed for transferring patient.

FAQ

1. which products you product?

Medical bed /medical trolley/home care bed / Electric Adjustable Beds

2. How do you ensure the products quality?



Production/inspection report before your order packed. we have obtained ISO Certificate.

3. How about your after-sale service?

We have a professional team to responsible for your after-sale service. Once the order confirmed, we have professional service follow this order.

4. What is your minium order quantity? If only buy one set to check, whether you can offer the best discount price for me?

Our Minium order quantity is one set. we will offer the best competitive price for each customer.

What is your delivery time? if the delivery time is urgent, how do you deal with?

Within 15 working days after order confirmed. If your delivery time is very urgent, we can delivery goods in 7 days. We take customer satisfying as our first goal.

5. What's the payment term?

Usually we accept T/T, L/C, Paypal, and Wester Union.

Warranty: 1-2 years due to different products, and lifetime maintaining.

A. Offer door to door service;

B. We send you free parts by courier and technical guidance by video/call;

C. We support you free training.