



Three-Function Electric Hospital Patient Bed

Welcome to wholesale and buy customized Three-Function Electric Hospital Patient Bed from our factory. You can rest assured to buy discount product from us. We look forward to cooperating with you, if you want to know more, you can consult us now, we will reply to you in time!

Product Description

Three-Function Electric Hospital Patient Bed

Three-Function Electric Hospital Patient Bed Made in China can be bought with Low Price from Shuaner, which is a professional high quality Products Manufacturers and Factory in China. We are a China company that can provide Customized Products service. If you want Pricelist and Quotation, you can ask us by leaving a message. We will get back to you as soon as possible.


High Grade 4pcs ABS Handrail ICU Household Nursing Used Three Functions Electric Bed For Hospital sale

1. PP headboard&foot board, 4pcs PP guardrail.
2. 4pcs 24V DCElectro-motor, 1pc control box.1pc handset
3. 4pcs urine hooks, 1pc IV pole,4pcs luxury castor (2pcs with brake ,2pcs without brake)
4. bed is special design for old patient caring, electric control it like Chair

Dedicated to strict quality manage and thoughtful buyer company, our experienced workers members are generally available to discuss your specifications and make certain full consumer gratification for Quality Inspection for China Three-Function Electric Hospital Patient Bed Manufacturers & Suppliers, "Change with the much better!" is our slogan, which means "A greater globe is before us, so let's take pleasure in it!" Change for that better! Are you completely ready?

Quality Inspection for China Three-Function Electric Hospital Patient Bed Manufacturers & Suppliers, Hospital Patient Bed, We have established long-term, stable and good business relationships with many manufacturers and wholesalers around the world. Currently, we are looking forward to even greater cooperation with overseas customers based on mutual benefits. Please feel free to contact us for more details.

Product Display



Packaging &Details

Packaging	20GP: 30 PCS 40HQ: 50 PCS
Warranty	2 Years
Delivery time	20 Days after Receiving 30% T/T In Advance

Product details

Model	Da-3a
Overall Size	2190*990*450-730mm
Bed Frame	Treated By Electro-Coating And Powder-Coating
Headboard / Footboard	PP Head / Footboard
Bedboards	Cold-Rolled Steel Bedboards
Handrails	PP Handrails
Hand Controller	Remote
Motor	Quiet And Robust Rlectric Actuators Provide Reliable Operation
	Remote Control Of all Electrically Operated Functions And Movements
	Linak Motor From Denmark And Jiechang From China
Bed Base	Steel Frame
Wheels	4 Double - Side Silcent wheels (150mm)With central locking system
Load Bearing	Fully Tested Sturdy Construction Capable Of Taking A Maximum User Weight Of Up To 300kg

Function

Backrest Max Upward Angle	75°
Footrest Max Upward Angle	40°



Height Adjustment	450- 730 mm
-------------------	-------------

Optional Accessories

I.V pole, Nurse control, CPR, Battery, Mattress

FAQ

1. which products you product?

Medical bed /medical trolley/home care bed / Electric Adjustable Beds

2. How do you ensure the products quality?

Production/inspection report before your order packed. we have obtained ISO Certificate.

3. How about your after-sale service?

We have a professional team to responsible for your after-sale service. Once the order confirmed, we have professional service follow this order.

4. What is your minium order quantity? If only buy one set to check, whether you can offer the best discount price for me?

Our Minium order quantity is one set. we will offer the best competitive price for each customer.

What is your delivery time? if the delivery time is urgent, how do you deal with?

Within 15 working days after order confirmed. If your delivery time is very urgent, we can delivery goods in 7 days. We take customer satisfying as our first goal.

5. What's the payment term?

Usually we accept T/T, L/C, Paypal, and Wester Union.

Warranty: 1-2 years due to different products, and lifetime maintaining.

A. Offer door to door service;

B. We send you free parts by courier and technical guidance by video/call;

C. We support you free training.