



### Three Function Manual Medical Hospital Bed

Welcome to wholesale and buy customized Three Function Manual Medical Hospital Bed from our factory. You can rest assured to buy discount product from us. We look forward to cooperating with you, if you want to know more, you can consult us now, we will reply to you in time!

#### Product Description

Three Function Manual Medical Hospital Bed made in China can be bought with low price from shuaner. It is a professional high quality Medical Bed manufacturers and factory in China. We are a China company that can provide customized products service. If you want pricelist and quotation, you can ask us by leaving a message. We will get back to you as soon as possible.

Three function manual medical hospital bed crank bed cheap beds

Advanced medical motor system(4 pcs motors,1 pcs control box, 1 pcs handset)

Urine hooks,I.V drip

Bed board with flexible connection

Central-controlled castor brake system

That has a sound business credit history, outstanding after-sales service and modern producing facilities, we have earned an superb popularity amid our buyers across the planet for Cheapest Price China Three Function Manual Medical Hospital Bed Manufacturers & Suppliers, "Change for the better!" is our slogan, which means "A better world is before us, so let's enjoy it!" Change for the better! Are you ready?

Cheapest Price China Three Function Manual Medical Hospital Bed, Medical Instrument, We supply experienced service, prompt reply, timely delivery, excellent quality and best price to our customers. Satisfaction and good credit to every customer is our priority. We focus on every detail of order processing for customers till they have received safe and sound items with good logistics service and economical cost. Depending on this, our products and solutions are sold very well in the countries in Africa, the Mid-East and Southeast Asia. Adhering to the business philosophy of 'customer first, forge ahead', we sincerely welcome clients from at home and abroad to cooperate with us.

#### Product Display



### Features

1. Detachable bedside and tailstock. With lock catach type system,it is easy to detach for rescue,special nursing,move the patients safty.
2. Handheld linetype remote control with large icon can freely handle.;groove on the guardrail can hold remote control to avoid to slide
3. The bed surface is shaped one for all by cold rolled steel plate.electrophoresis and powder flame plating ensure inner of the pipe is protected by paint and the outer and inner quality is same
4. Imported chained motor ,VO level firproof material is very smooth ,it is easy for clean.the motor's noise is very low and lifetime is long.
5. Back side and leg side is connected,when the back side up and down,the leg side will coordinate,thus can avoid bed sore and reduce medical worker's work pressure.

### Packaging&Details

Packaging	20GP: 50PCS 40HQ: 80PCS
Warranty	2 Years
Delivery time	20 Days after Receiving 30% T/T In Advance

### Product Details

Model	SAE-YB-3a
Overall Size	2190*980*400-630mm
Bed Frame	Treated By Electro-Coating And Powder-Coating
Headboard / Footboard	PP head&foot board
Bedboards	Cold-Rolled Steel Bedboards
Handrails	PP guardrail
Hand Controller	Cranks
Bed Base	Steel Frame
Wheels	4 Double - Side Silcent wheels (150mm)
Load Bearing	Fully Tested Sturdy Construction Capable Of Taking A Maximum User Weight Of



	Up To 300kg
--	-------------

### Function

Backrest Max Upward Angle	75°
Footrest Max Upward Angle	40°
Height Adjustment	400 - 630 mm

### FAQ

1. which products you product?

Medical bed /medical trolley/home care bed / Electric Adjustable Beds

2. How do you ensure the products quality?

Production/inspection report before your order packed. we have obtained ISO Certificate.

3. How about your after-sale service?

We have a professional team to responsible for your after-sale service. Once the order confirmed, we have professional service follow this order.

4. What is your minium order quantity? If only buy one set to check, whether you can offer the best discount price for me?

Our Minium order quantity is one set. we will offer the best competitive price for each customer.

What is your delivery time? if the delivery time is urgent, how do you deal with?

Within 15 working days after order confirmed. If your delivery time is very urgent, we can delivery goods in 7 days. We take customer satisfying as our first goal.

5. What's the payment term?

Usually we accept T/T, L/C, Paypal, and Wester Union.

Warranty: 1-2 years due to different products, and lifetime maintaining.

A. Offer door to door service;

B. We send you free parts by courier and technical guidance by video/call;

C. We support you free training.